

HOW TO USE THE COMPRESSOR-DRIVEN NEBULIZER

Assembling the Nebulizer:

1. Wash your hands.
2. Remove all nebulizer parts from the plastic bag, including mask/mouthpiece, flex tubing, nebulizer cup, and connecting tubing.
3. Assemble the nebulizer by attaching the t-piece to the nebulizer cup, and then attaching flex tubing to one end of the t-piece and mouthpiece or mask to the other end.
4. Attach connecting tubing to the outlet on the bottom of the nebulizer cup. Attach the other end to the outlet on the nebulizer compressor.

Using the Nebulizer:

1. Sit upright in a comfortable chair. An infant or small child may be held upright on lap or in arms.
2. Place mouthpiece in mouth or place mask on infant or child.
3. Turn switch on compressor to the “on” position.
4. Breathe normally through the mouthpiece or mask, taking a deep breath every minute or so. Continue until neb cup begins to sputter. Shake or tap the cup and continue breathing through the mouthpiece/mask until nearly all liquid in neb cup is gone.

